



Description

For use as an emergency call button

Theory of Operation

The Emergency Button provides patients with a means to call care providers for assistance. The large spring-loaded button acts as a momentary contact switch (normally open or normally closed depending on model). Once activated, the button transmits a signal to the owner's monitoring system and provide the patient with an optional audio message. The ON/OFF switch on the underside of the unit can enable and disable the audio speaker. To power the system off, remove the batteries from the underside of the unit.

Additional Features

A wireless transmitter may be incorporated. A switch on the base of the unit can disable the speaker without disabling the wireless signal.

Construction

Parameter	Typical Val	Unit
Housing Diameter	3.6	in
Housing Height	1.45	in
Max Operating Current	500	mA
Max Operating Voltage	3.3	V

Specifications

Parameter	Description
Interface	Normally open and normally closed momentary contact switches available.
Audible Message	Can be customized

Physical Dimensions

Side View

